

# Perfecting Performance



Occlutech® **Sizing Balloon**




## Product Range - Dimensions and sizes

Art No.	Balloon Ø	Fill volume	Max. fill volume	Ø Balloon at max. fill volume	Balloon length
53SB254	25 mm	25 ml	35 ml	27 mm	45 mm
53SB355	35 mm	60 ml	90 ml	40 mm	50 mm

### State-of-the-art Sizing

- Rapid inflation and deflation
- Soft compliant balloon
- Atraumatic tip



#### Designations:

- (1) Usable length
- (2) X-ray markers
- (3) Soft Tip
- (4) Inflation port

Art No.	Balloon Ø	Balloon length	Catheter length	Diameter	Guide wire
53SB254	25 mm	45 mm	700 mm	8F	0.035" (Inch)
53SB355	35 mm	50 mm	700 mm	8F	0.035" (Inch)

#### Customer Service International

Occlutech, Rönnowsgatan 8, S-252 25 Helsingborg, Sweden  
 Tel +46 42 33 65 25, Fax +46 42 311 09 70  
 order@occlutech.com, [www.occlutech.com](http://www.occlutech.com)

