


**Perfecting** Performance



**Occlutech® Stiff Guide Wire** The Obvious Leader



# Occlutech® Stiff Guide Wire

## The Obvious Leader

The Occlutech Stiff Guide Wire has been developed using a stainless steel core and a unique PTFE coating. This combination gives well-balanced support for day-to-day cases, as well as advanced ones.

A soft, atraumatic J-tip, in combination with the exceptional torque control, allows exceptional navigation in difficult and challenging geometry.

This stiff configuration coupled with the flexible tip allows a good push while making the wire very safe and easy to use.

- **Soft atraumatic J-tip** – allowing safe navigation
- **Stiff configuration** – supporting easy entry into the left atrium
- **Completely coated** – minimizing the risk of thrombus in the left atrium



## Ordering information

### Occlutech® Stiff Guide Wire

Article No.	Product Name	Ø (inch)	Length (cm)	Tip Shape (mm)
52GWR01	Occlutech Stiff Guide Wire	0.035	260	J-3mm

**Customer Service International**  
Occlutech, La Cours Gata 2, S-252 31 Helsingborg, Sweden  
Tel +46 (0)42 33 65 32, Fax +46 (0)42 311 09 70  
order@occlutech.com, [www.occlutech.com](http://www.occlutech.com)

© Occlutech 2015. All rights reserved. P07F06.020.03

