

Perfecting Performance



Occlutech Sizing Balloon



OSB

Occlutech Sizing Balloon

OSB

Changing the dimension

The Occlutech Sizing Balloon is intended for venous access. It is used without a guide catheter and is introduced directly via the guide wire.



Designations:

- (1) Usable length
- (2) X-ray markers
- (3) Soft tip
- (4) Inflation port



Product specifications

| Article No. | Balloon Ø Nominal [mm] | Nominal Fill Vol. [ml] | Max. Fill Vol. [ml] | Balloon Ø Max. Fill Vol. [mm] | Balloon Length [mm] | Usable Length [mm] | Guide Wire [inch] |
|-------------|------------------------|------------------------|---------------------|-------------------------------|---------------------|--------------------|-------------------|
| 53SB254 | 25 | 25 | 35 | 27 | 45 | 700 | 0.035 |
| 53SB355 | 35 | 60 | 90 | 40 | 50 | 700 | 0.035 |



Occlutech International AB
Landskronavägen 2, SE-252 32 Helsingborg, Sweden
order@occlutech.com, www.occlutech.com

